

Inspiring innovation.

# ASSIMILABLE AMINO NITROGEN PROCEDURE RAPID METHOD 1998

### Equipment

- Spectrophotometer, UV at 335nm, zero on DI water
- Cuvettes, 10mm
- Pipettor, 50μL
- Adjustable micro-pipettor with disposable syringes
- 50 mL centrifuge tubes
- Centrifuge

#### **Reagents**

- NAC Buffer (Reagent buffer-No OPA)
- OPA Solution, 5%
- 10mM Isoleucine Standard (10mM IIe)
- DI Water

### Procedure

Mix reagents as needed:

**NAC Buffer**, pour premeasured crystalline NAC into premeasured borate solution. Use within 2 days of mixing. Refrigerate for stability and use at room temperature.

5% OPA, add the contents of ethanol vial (5 mL) to the premeasured OPA. Use within 1 day.

#### STANDARD CURVE

1. Add 3 mL NAC buffer into cuvettes. Add 10mM lle standard and water as follows:

Cuvette	1	2	3	4
10mM lle	0	10µL	30µL	50µL
DI Water	50µL	40µL	20µL	0

- 2. Add 50 $\mu\text{L}$  5% OPA solution to each cuvette.
- 3. Read Absorbance values and record in log book.

The indications supplied are based on our current knowledge and experience, but do not relieve the user from adopting the necessary safety precautions or from the responsibility of using the product(s) properly.



Inspiring innovation.

## Samples

Setting up calculations in EXCEL is most expedient. Note: To convert to ppm Nitrogen, the multiplication factor is 14.

1. Key in (on HP scientific calculator)

f	GSB	(clears statistical register)	
Concentration value 1 ( y ), Enter			
Absorbance value 1 ( x ) $\Sigma$ +			
Concentration value 2, Enter			
Absorbance value 2 $\Sigma$ +			
Etc		(Enter all the concentration (y) and absorbance (x) values)	
$f E \Sigma E$		(Calculates the linear regression equation and displays the y-intercept)	
$X \Leftrightarrow Y$		(Displays the slope)	
f∙,Ethen X	⇔Y	(Displays r value)	
Sample Absorbanc	te value $f \bullet$	(Displays concentration value) (multiply by 14 to convert to ppm Nitrogen)	

## Disposal

Dispose in sink.

The indications supplied are based on our current knowledge and experience, but do not relieve the user from adopting the necessary safety precautions or from the responsibility of using the product(s) properly.

Revision: August 2019