





Unique Features

- Electronically regulated heating element—you can walk away and the boiling point remains steady.
- **Pyrex boiling chamber**—allows you to see the sample and it's easy to clean.
- Circulating condenser—protects the sample from overheating or boiling away.
- 0.1% vol/vol precision—determination of alcohol content in dry wines.
- Faster analysis time—each sample takes less than 6 minutes after the initial determination for the boiling point of water. With continuous cooling water flow through the condenser and easy sample introduction and waste disposal, analysis time is reduced.
- **Ready to use**—requires electrical hook up, access to cooling water for the condenser and a drain for sample disposal. Each unit comes complete with certified thermometer and a sliding ruler for alcohol calculation.

Description and Use

The electric ebulliometer has a built-in heating regulator. During the first 3 minutes, the heating element is on full power. After 3 minutes, the power is cut by half to allow the boiling point to stabilize. After stabilization, the boiling point will be maintained for an extended period of time.

The indications supplied are based on our current knowledge and experience, but do not relieve the user from adopting the necessary safety precautions or from the responsibility of using the product(s) properly.

Revision: September 2019