

Inspiring innovation.

OBJECTIVE MICROSCANS

Bacteria:	Average count/views
Heavily loaded	120 or more cells/view 20 or more groups/view
Very large population	36-120 cells/view 10-20 groups/view
Large population	11-35 cells/view 5-10 groups/view
Moderate population	4-10 cells/view 1-4 groups/view
Small population	1-3 cells/view 5-30 groups/30 views
Few	<1 cells/view <5 groups/30 views

Yeast:	Average count/views
Heavily loaded	over 50 cells/view
Very large population	21-50 cells/view
Large population	8-20 cells/view
Moderate population	3-7 cells/view
Small population	1-2 cells/view
Few	<1 cell/view

The indications supplied are based on our current knowledge and experience, but do not relieve the user from adopting the necessary safety precautions or from the responsibility of using the product(s) properly.

Revision: March 2020