

TARTARCHECK PLUS

RAPID COLD STABILITY MEASUREMENT

Quick Facts

What: Determination of Cold Stability (Tartrate Stability) of wines.

How: Measures conductivity drop after contact with KHT.

Why: Allows for customers to determine readiness for bottling and validation of any cold treatments.



What's Included

- Tartarcheck Plus (no computer necessary)
- Sample Syringe
- KHT Syringe
- KHT

Further Considerations

- No computer requirement (which means no software updates)
- Uses Δ (change) in conductivity to predict stability

Information Resources

- [Enartis Product Page](#)
- [Manufacturer Website](#)
- [Operational Video](#)

Purchasing Information

Product	Item Number
Tartarcheck Plus	50-700-1000
KHT	10-207-0100 or 10-207-0500
Syringe	20-219-0060